

Slide Dryer TSD-4E

 Thermoline


The Thermoline Slide Dyer is a High Capacity Slide Dyer with Advanced Controls. The extra-large capacity and angled slide storage area will make it an essential part of your histology workflow. The easy to use touch screen interface allows for accurate temperature control for ideal slide preparation.



Features

- 56 slide capacity (conventional slides): 7 blocks and 8 slides on each block.
- High-definition touch screen and user-friendly interface for intuitive operation.
- Microprocessor controlled device with precise temperature performance that can keep the temperature accuracy within 1°C.
- Advanced safety over-temperature protection that shuts off the devices when the temperatures is over 70°C.
- A nonvolatile memory retains the last temperature set point even after powering off.
- Timer and alarm reminder.
- Auto turn-on/ turn-off function eliminates wait time.

Slide Dryer



Specifications

Product Name	Slide Dryer
Model	TSD-4E
Temperature Setting Range	0°C~99°C (1°C increment)
Voltage	AC230V±10%, 50Hz±1Hz
Power	200W
Dimensions (L x W x H)	327mm x 340mm x 166mm
Maximum Slide Capacity	56
Temperature Range	Ambient -99°C
Temperature Accuracy	±1°C
Over-temperature Alarm	YES
Surface Size	34 x 22cm
Display	LED Display
Weight	8.5kg

